


Antivirové šátky nanoSPACE

nanoSPACE®

Kategorie zařazení: Komunitní nano rouška

Výrobce nanoSPACE s.r.o. prohlašuje, že Antivirový šátek nanoSPACE zakoupený na www.nanoSPACE.cz splňuje následující parametry:

- kategorie produktu: Komunitní nano rouška
- použitý nanovláknenný materiál: 95 % polypropylen a 5 % PA6 se zachytem submikronových částic 97-99,9 %
- původ materiálu: Pardam Nano4fibers, Roudnice nad Labem
- výsledky nezávislého měření dTest (12/2020): filtrační účinnost FFP2
- test zachytu virů (VFE) v Nelson Labs v USA: zachyt více než 99,9 %
- Výrobce je hrdým členem Asociace nanotechnologického průmyslu ČR a je součástí projektu Česko je nano



Results:

Test Article Number	Percent VFE (%)
1	>99.9 ^a
2	>99.9 ^a

^a There were no detected plaques on any of the Andersen sampler plates for this test article.

The filtration efficiency percentages were calculated using the following equation:

$$\% VFE = \frac{C - T}{C} \times 100$$

C = Positive control average

T = Plate count total recovered downstream of the test article
Note: The plate count total is available upon request

Lucie Konečná
provozní ředitelka nanoSPACE