nano SPACE®

Rohova 98, Dolejší Předměstí, 344 01 Domažlice 377 311 311 info@nanospace.cz www.nanospace.cz

Washable Designer Nano Face Masks by nanoSPACE

Classification: Community nano face mask

nanoSPACE s.r.o. as the manufacturer hereby declares that the nanoSPACE Washable Designer Nano Face Mask purchased on www.nanoSPACE.cz complies with the below parameters:

- Nanofiber material: 95% polypropylene, 5% PA6 with submicron particle filtrationefficiency of 97-99.9%
- Material supplied by: Pardam Nano4fibers, Roudnice nad Labem, Czech Republic
- Virus filtration efficiency (VFE) tested by Nelson Labs, USA: test proven filtrationefficiency exceeding 99.9%

The manufacturer is a proud member of the **Czech Nanotechnology**

Industries Association and a partner in the Czech Is Nano project

Results:

Test Article Number	Percent VFE (%)
1	99.5
2	99.0
3	>99.9ª

° There were no detected plaques on any of the Andersen sampler plates for this test article. The filtration efficiency percentages were calculated using the following equation:

% VFE =
$$\frac{C-T}{C} \times 100$$

versited nenoSPAC workhoble designer mene mask with the second C = Positive control average

 $\mathsf{T}=\mathsf{P}\mathsf{late}$ count total recovered downstream of the test article Note: The plate count total is available upon request

Lucie Konečná ditelka nanoSPACE provoz